



PRODUCT CATALOGUE

CERTIFIED QUALITY

CIZETA MEDICALI S.p.A. products comply with the main international quality standards for the protection and safeguarding of patients. Product quality is the key to its success on a national and international scale.



UNI EN ISO 9001 e UNI CEI EN ISO 13485

CIZETA MEDICALI has a quality management system compliant with the UNI EN ISO 9001 and UNI CEI EN ISO 13485 standards. Through these certifications, obtained in 2001, CIZETA MEDICALI guarantees the highest quality standards of its products and service.



MDR 2017/745
Class I Medical Device

CE Mark

CIZETA MEDICALI products are Class I Medical Devices and satisfy the requirements of MDR 2017/745.



OEKO-TEX STANDARD 100

All products manufactured by the Group are OEKO-TEX certified, to guarantee the absence of harmful substances to safeguard Patients' health.



RAL

Certification issued by the Hohenstein Research Institute (Germany), which guarantees the quality of the stockings through periodic quality testing in compliance with the RAL-GZ 387 standard.



ASQUAL

Certification issued by the ASQUAL Quality Association (France), which guarantees that the certified products comply with specifications of the NF G30-102B standard, guaranteeing their quality.



CFDA

Certification issued by the CFDA (Chinese Food and Drug Administration) (China), which guarantees the conformity of the certified products with the specifications of the YY/T 0853-2011 standards.



FDA

The Food and Drug Administration (FDA) is an agency of the United States Department of Health and Human Services. CIZETA MEDICALI products are approved by the FDA, whose mission is to protect public health, safety and efficacy.



EUROPEAN MANUFACTURERS
FEDERATION FOR COMPRESSION
THERAPY AND ORTHOPAEDIC
DEVICES

CIZETA MEDICALI is a founding member of Eurocom, the European association of manufacturers for compression therapy and orthopaedic devices. Only companies that produce medical compression stockings with RAL-GZ 387 certification can be members of Eurocom. The mission of the association is to promote awareness of the quality of RAL-GZ 387 certified stockings among industry operators and patients.

THE COMPRESSION CLASSES

CIZETA MEDICALI S.p.A. specialises in the manufacture of elastic compression braces used to treat phlebological and lymphatic diseases.

Elastic compression certainly remains the best known treatment - effective and easily adapted to such diseases.

The application of external compression can at least partly counterbalance the effects of haemostasis and lymphostasis.

The most suitable compression class is chosen according to the stage of the disease.



The braces are produced in four compression classes, expressed in mmHg and based on the compression exercised at the ankle.



Product certification according to the German RAL-GZ 387 standards



Product certification according to the French NF G30-102B standards

INTERNATIONAL PRODUCT CERTIFICATIONS		
Compression class		
Ccl. 1	18-21 mmHg	10-15 mmHg
Ccl. 2	23-32 mmHg	15-20 mmHg
Ccl. 3	34-46 mmHg	20-36 mmHg
Ccl. 4	>49 mmHg	>36 mmHg

The most well-known international certification bodies are the Hohenstein Research Institute (Germany), which guarantees the quality of the stockings through periodic quality testing of compliance with RAL-GZ 387 regulations, and the ASQUAL Quality Association (France), which guarantees that the certified products comply with NF G30-102B regulations.

The benefits of certified medical compression stockings are:

- correct compression guaranteed at the ankle, in accordance with the declared compression class
- correct graduated compression profiling (gradual decrease of the pressure exerted from the ankle upwards)
- six-month guarantee of compression and graduation, from their first use



Characterized by excessive accumulation of the adipose tissue, Lipedema is a chronic inflammatory disorder affecting predominantly the feminine sex.

According to its appearance and location area, generally excluding the extremities of the body, it can be classified into stages and types.

Typical symptoms of this condition are oedema, heaviness and pain to the touch that, if neglected, can lead to a degeneration of Lipedema towards the higher stages, causing walking difficulties.

LIPEDEMA LINE



Its treatment requires the daily use of a elastocompressive garment, able to counteract oedema and reduce painful symptoms, accompanied by an healthy lifestyle (diet and physical activity).

LIPOCARE[®] LIGHT

Lipedema line
Circular knit medical compression garment




CIRCULAR KNIT



STANDARD GARMENTS

CIRCUMFERENCE (cm)							
Size	1/XS	2/S	3/M	4/L	5/XL	6/XXL	
Thigh	G	45-55	50-60	55-65	60-71	66-78	73-85
Calf	C	29-35	31-38	34-41	37-44	40-48	44-52

LENGTH (cm)		
	Short	Normal
AG	66-73	74-82



COMPRESSION CLASS		
REF.	23 - 32 mmHg	ARTICLE
	L4ECT	Short panty under the knees (CT)
	L4NCT	Normal panty under the knees (CT)
	L4EBT	Short leggings (BT)
	L4NBT	Normal leggings (BT)



LIPOCARE® LIGHT

Lipedema line

Circular knit medical compression garment



LIPOCARE® is the line of compression garments dedicated to the treatment of lipedema of lower limbs.

The LIPOCARE® LIGHT line, manufactured with medical circular knit technology and high-quality yarns, meets the therapeutic needs of early-stage lipedema by reducing the feeling of heaviness and pain in legs, and ensuring a feeling of wellness.

The special rib pattern has a draining effect and bestows the product a modern and fashionable look, ideal for everyday wear.

Available as leggings (BT) and panty under the knees (CT).

CHARACTERISTICS

- The medical design allows the correct preservation of compression levels ideal to relieve the symptoms of lipedema
- The high-quality microfibre yarn makes the product extremely soft to the touch and gentle on the skin, and ensures high levels of breathability
- With knitted-in strap that allows correct positioning of the product
- Perfect fit, form-fitting
- Easy to put on

INDICATIONS

CLASS II (23 - 32 mmHg)

Lipedema, lipolymphedema, primary and secondary lymphedema in the various clinical developing stages, edematous-fibro-sclerotic-panniculopathy (known as "orange peel skin").

CONTRAINDICATIONS

Acute dermo-hypoderma, acute skin infections, obliterative peripheral arteriopathy, known contact allergy to elastic braces.

COMPOSITION

Polyamide 60%, elastane 40%.



WIDE RANGE



BREATHABILITY

COMPRESSION

Ccl. 2 23-32 mmHg



STANDARD
100
074156.0
CENTROCOT

COLOURS

862 black



LIPOCARE® STRONG

Lipedema line
Flat knit medical compression garment




FLAT KNIT



STANDARD GARMENTS

CIRCUMFERENCE (cm)							
Size	1/XS	2/S	3/M	4/L	5/XL	6/XXL	7/XXXL
Measurement points							
cT	65-75	75-85	85-95	95-105	105-115	115-125	125-135
cH	85-95	95-105	105-115	115-125	125-135	135-145	145-155
cG	44-55	48-59	52-63	56-67	60-71	64-75	68-79
cF	43-52	47-56	51-60	55-64	59-68	63-72	67-76
cE	32-37	35-40	38-43	41-46	44-59	47-52	50-55
cD	28-33	30-35	32-37	34-39	36-41	38-43	40-45
cC	30-36	33-39	37-43	40-46	43-49	46-52	49-55
cBⁱ	23-29	26-31	29-35	32-38	35-41	38-44	41-46
cB	18-21	20-23	22-25	24-26	25-29	28-32	31-35



LENGTH (cm)		
	Short	Normal
AG	66-73	74-82

COMPRESSION CLASS			
REF.	23 - 32 mmHg	ARTICLE	COLOURS
	FM4ECT	Short panty under the knees (CT)	604 862
	FM4NCT	Normal panty under the knees (CT)	
	FM4EBT	Short leggings (BT)	604 862
	FM4NBT	Normal leggings (BT)	



LIPOCARE® STRONG

Lipedema line

Flat knit medical compression garment



LIPOCARE® is the line of compression garments dedicated to the treatment of lipedema of lower limbs.

The LIPOCARE® STRONG line, manufactured with medical technology and high-quality yarns, meets effectively the most complex therapeutic needs and, thanks to its flat-knitted structure, has high containment capacities of adipose accumulation, helping to improve the patient's mobility.

Available as leggings (BT) and panty under the knees (CT).

CHARACTERISTICS

- The medical flat-knitted design ensures high levels of containment and promotes mobility
- The high-quality microfibre yarn makes the product extremely soft to the touch and gentle on the skin, and ensures high levels of breathability
- The weave in the groin area is lightened to facilitate the correct positioning of the product
- The reduced levels of compression at the ankle avoid the tourniquet effect
- With adjustable velcro waist strap for optimal fit and maximum support

INDICATIONS

CLASS II (23 - 32 mmHg)

Lipedema, lipolymphedema, primary and secondary lymphedema in the various clinical developing stages, phlebolympedema, edematous-fibrosclerotic-panniculopathy (known as "orange peel skin").

CONTRAINDICATIONS

Acute dermo-hypoderma, acute skin infections, obliterative peripheral arteriopathy, known contact allergy to elastic braces.

COMPOSITION

Polyamide 74%, elastane 26%.



WIDE RANGE



BREATHABILITY

COMPRESSION

Cl. 2 23-32 mmHg



COLOURS

604 havana



862 black





SCIENCE & MEDICAL
QUALITY

CIZETA MEDICALI S.p.A.
Via IV Novembre, 46 - 20012 Cuggiono (MI) Italy
Tel. +39 02.9721811 - Fax + 39 02.97240734
www.cizetamedicali.com - info@cizetamedicali.com



@cizetamedicali